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Inertsil® HPLC COLUMNS

Operating Instructions

Thank you for purchasing the Inertsil HPLC column.

Inertsil packing materials are subject to a rigorous array of QC tests in the ISO9001 compliant facility, with special emphasis on reagent purity, raw material traceability and consistency in raw materials, and final products. A detailed analysis of all of the physical and chemical properties of the Inertsil column, combined with tests for chromatographic selectivity and column packing material efficiency, ensure that each lot of Inertsil columns is identical to all previous lots and column-to-column reproducibility is of the highest order.

To maintain and maximize peak performance of Inertsil HPLC column, and ensure long life and stability of columns please read the following instructions before use.

(1) UNPACKING:

☑ Check if anything is missing or damaged. If there are any signs of damage, notify your local GL representative at once.

(2) SHIPPING SOLVENTS:

☑ Inertsil columns are shipped with the solvent used for the final QC test of the column, as detailed on the test chromatogram delivered with the column. When switching between solvents with vastly different polarities, it may be necessary to first purge the column with a mutually miscible solvent such as Isopropyl alcohol or Dioxane at a reduced flow rate (approximately 50% of normal). Pump that solvent for five times more column volume (e.g. 10mL for a 150 x 4.6mm I.D. column).

Type	Principal columns	Shipping Solvent
Normal-phase	Inertsil SIL-100A, CN-3, NH2, Diol...etc.	n-Hexane : Ethanol
Reversed-phase	Inertsil ODS-4, ODS-3, C8-4,C8-3, Ph-3...etc.	Water : Acetonitrile Water : Methanol
HILIC-phase	Inertsil HILIC, Amide	Water : Acetonitrile
Ion-exchange	Inertsil CX, AX...etc.	Methanol 100%

Note: With Inertsil SIL-100A, it is time consuming to flush the buffer or water in the pores of the packing after using reversed-phase eluent. If a subsequent analysis in the normal phase mode is performed, and the flush has not been thoroughly completed, the elution order and retention time may vary significantly. Therefore this type of column should be used solely for one type of analysis.

(3) HANDLING:

☑ Shock : Do not drop or bump columns, to avoid a deterioration of the column performance.

☑ Eluent pH : Inertsil columns provide optimum lifetime when used with eluents within the pH range of 2.0 to 7.5. Higher pHs leads to dissolve the silica gel and lower pHs cause to hydrolyze away some of the bonded phase. While many customers use the columns outside both sides of the pH range with excellent results and good column lifetime, the best lifetimes are usually obtained at intermediate pH conditions.

☑ Pressure : To maximize column life, use the columns within the pressure shown in the following table.

Particle Size	Recommended Operating Pressure (Maximum Operating Pressure)
2 µm, UHPLC PEEK2 µm	80 (MPa)
3 µm HP, UHPLC PEEK3 µm	50 (MPa)
3 µm, 5 µm, PEEK5 µm	20 (MPa)

☑ Avoid rapid pressure fluctuation.: Make sure to disconnect the column from the system after confirming the display of the pressure gauge showing zero "0" value.

☑ Sample Dissolution : Samples should be dissolved in the eluent or solvent weaker than the eluent, which helps avoid sample precipitation at the column inlet/head and inconsistent retention values.

☑ If peaks are tailing more on the early eluting compounds than later eluting compounds, there is a possibility that there is a dead volume. In such case, check all column connections are properly connected.

Also, use appropriate internal diameter and length size of tubing at the injector and detector, especially when using semi-micro size columns to avoid system dead volumes.

☑ The most common cause of columns back pressure increase or double peaks is blockage of the inlet filter by sample particulates, particles created by aging pump seals, or large quantities of lipophilic compounds adsorbing to the head of the column.

• Pass the eluent and samples through a GL-Chromatodisk or other quality 0.45µm membrane filter before use.

• Installation of a GL Sciences' Preclean ORG between the pump and injector is also effective to remove particulates or contaminants in the mobile phase.

• Make sure to filtrate the sample using "GL Chromatodisc" (syringe filter) before injecting to the column.

• The HPLC column can be also be protected from contamination by using "Guard Column for UHPLC" or

"Cartridge Guard Column E", which results in longer column life and more consistent analytical results.

☑ When observing baseline noise, make sure HPLC-grade solvents are used. Aged detector lamp could be another possible cause.

Make sure the solvents are fully degassed and free from air bubble to prevent bubble formation in the pump and detector.

☑ One of the most common sources of poor performance is system dead volume, particularly between the injector and the column inlet.

Be certain that dead volume is minimized by using proper fittings and tubes.

☑ Baseline drift and noise can be caused by defective pumping eluent due to air bubbles or reduction of light intensity when using a UV detector. Bubbles can form in the detector flow cell if the eluent is not degassed properly before introduction into the column.

About Inertsil Amide columns

☑ Always maintain at least 50 % organic solvent (e.g.: acetonitrile) in your mobile phase or gradient. Operating with aqueous concentrations greater than 50 % will result in decreased column lifetime.

☑ To maximize performance under HILIC mode, use acetonitrile/ammonium acetate aqueous solution = 80/20. Stronger retention can be obtained by increasing the concentration of organic solvent, which is the feature of a HILIC mode. (e.g.: acetonitrile/ammonium acetate aqueous solution = 90/10)

About PEEK (poly ether ether ketone) columns

☑ End-fittings are 1/16" type. (No. 10-32 UNF, Tubing Tip Length 2.4 mm)

☑ Use PEEK fittings when connecting columns. Do not overtighten since this may break the screw and lead to column breakage.

☑ Avoid the use of THF (tetrahydrofuran) and chloroform as it can weaken the PEEK hardware and cause it to become brittle.

(4) STORING OF COLUMNS:

• After using the columns with eluent containing buffer or ion-pair reagent, wash the column thoroughly with a salt-free eluent before storing.

• When storing the column for a long period, store it with 100% organic solvent.

• Seal the column with the plugs provided and store it with a stable temperature.

Inertsil columns are manufactured, inspected, packaged and shipped under strict standards of quality control. Should you find any defect in performance, please contact us or your GL Sciences' provider, who will ensure your complete satisfaction. We regret that we cannot guarantee the lifetime of columns, nor can we accept any claim when their performance has deteriorated due to no-compliance with the above operating instructions or as a result of normal aging.

For more information visit www.glsciencesinc.com and view detailed Column Care and Use information.

GL Sciences Inc.

22-1 Nishishinjuku 6-chome, Shinjuku-ku,
Tokyo 163-1130. Japan
Tel. +81-3-5323-6620 Fax. +81-3-5323-6621
<http://www.glsciences.com>
E-mail:world@gl.co.jp