

Easy Concentration & Analysis of Fresh and Spoiled Cabbage by MonoTrap[®] - Thermal Desorption

Using MonoTrap RGC18TD containing graphite carbon, the volatile compounds of fresh and spoiled cabbage are determined by screening analysis.

Protocol

Cabbage

Put 25g of chopped cabbage into a 100mL vial

Spoiling Stage

Leave it at 60 °C for 3 – 7 days

Passive Sampling
MonoTrap RGC18TD × 2pcs

Room temperature
For 3 hours

TD-GC/MS-O




<Left>
Freshly
chopped

<Right>
After 3 days at
60°C

GCO Sniffing Port (OP275)

Spoiling Stage & Volatile Compounds




GC Conditions

System	: GC/MS-Thermal Desorption (T-Dex II)
Column	: InertCap Pure-WAX 0.25 mm I.D. × 60 m, df = 0.25µm
Col.Temp.	: 40°C (5 min) - 6°C/min - 250°C
Carrier Gas	: He, 1mL/min (constant flow)
Desorb Temp.	: 200°C
Time	: 5 min
Flow	: 5mL/min
Split	: Splitless
Cryo Trapping	: -150°C
Injection Temp.	: 250°C
Detection	: MS Scan (28.5 - 600 m/z)

Split the capillary column outlet and simultaneous measurement by GC/MS and GC/O. The sensitivity of MS is 1/10 of ordinal detectors. Due to the graphite carbon effect, SN compounds can be effectively detected.

Organoleptic Data of Volatile Compounds from Spoiled Cabbage



Compound	Smell	Compound	Smell
1. Dimethyl disulfide	Spoiled egg	11. Acetic acid	Sour
2. Hexanal	Fresh	12. 1-Octen-3-ol e alcohol	
3. 2-Heptanone		13. Hexane, 1-nitro-	Metal
4. 1-Butanol, 2-methyl-	Raw garbage	14. Benzaldehyde	
5. 1-Butanol, 3-methyl-	Bitter	15. 4-Hexen-1-ol,	
6. 1-Pentanol	Fresh	16. Benzonitrile	Bitter
7. Hexanenitrile	Spoiled	17. S-Methyl methanethiosulphonate	Spoiled
8. 2-Butanone, 3-hydroxy-		18. (2,6,6-Trimethyl-2-hydroxycyclohexylidene)-acetic acid lactone	
9. Dimethyl trisulfide	Spoiled	19. Indole	Foul
10. 1-Hexanol	Irritating		