

Analysis of pesticide residues in agricultural products using automated Liner Exchange Difficult Matrix Introduction (LINEX-DMI) technique coupled to Gas Chromatography/ Mass Spectrometry (GC/MS)

Y. Takei^{a)}, K. Iijima^{b)}, M. Saka^{b)}



a) GL Sciences, Inc., 237-2 Sayamagahara, Iruma, Saitama, 358-0032, Japan E-mail: takai@gl.co.jp

b) The Institute of Environmental Toxicology (IET), 4321, Uchimoriya-machi, Joso-shi, Ibaraki, 303-0043, Japan



1. Introduction

For the determination of multiple of pesticide residues in foods, the complex preparation is generally required and, therefore, is time-consuming due to the complexity of both the matrices and the interfering compounds. A technique called the LINEX-DMI, in which automatic liner-exchanging (LINEX) is combined with difficult matrix introduction (DMI), has become a powerful tool for overcoming the problems. In the LINEX-DMI approach, the GC liner is automatically exchanged, and the non-volatile matrices such as pigments and lipids trapped in the liner can be removed at each run. With the application of the LINEX-DMI system, the sample clean-up procedures are simplified, and the contamination of the capillary column and the detector caused by non-volatile matrices is effectively eliminated.

A main goal of the study is to apply the LINEX-DMI system into the determination of the pesticides that are sensitive to thermal decomposition and adsorption.

2. LINEX-DMI technique

The combination of the automatic liner-changing with the DMI is very effective for the determination of the pesticides in foods and the PCB in insulation oil. When applying the LINEX-DMI into a capillary GC/MS system, the contamination caused by the direct introduction of dirty samples can be prevented.

- In the DMI system, a micro-vial is placed into the GC liner. (Figure 1).
- A liquid or solid sample is added into the micro-vial. (size: 2mmI.D x 40mm length ; Figure 1).
- The micro-vial is inserted and fixed on a narrow diameter of the liner. (Figure 1c).
- The liner inserted with the micro-vial is automatically installed into the LINEX injection port by the FOCUS autosampler. (Figure 2, 3, 4)
- A septum purge line can be performed using the LINEX head. (Figure 2) The PTV injection results in that the non-volatile matrices remain in the micro-vial, and prevents the contamination of the capillary column and MS. (Figure 1 e)




Figure 1. DMI technique




Figure 2. LINEX Head on OPTIC3 injector



Figure 3. LINEX Gripper on FOCUS robot



Figure 4. GC/MS system on LINEX system

3. Experimental

3.1. Chemicals and samples

- ◆ Internal standard component; Phenantrene - d10
- ◆ Standard sample 1; composite stock standard solution (1ng/μL) of 11 paraffinic hydrocarbons (n-C10, n-C12, n-C14, n-C16, n-C18, n-C20, n-C22, n-C24, n-C28, n-C30 and internal standard) was prepared in CHCl₃.
- ◆ Standard sample 2; composite stock standard solution (5ng/μL) of 15 pesticides (Fenobucarb, Bendaiocarb, Dimethoate, Lindane, Diazinon, Ethiofencarb, Carbaryl, Fenitrothion, Methiocarb, Aldrin, Fenthion, Dieldrin, Endrin, p,p'-DDT, Etofenprox and internal standard) was prepared in acetone.
- ◆ Standard sample 3; composite stock standard solution (1ng/μL) of 243 pesticides and internal standard was prepared in acetone.
- ◆ Extraction procedures are shown in Figure 10.

3.2. Instrumentation

3.2.1. LINEX-GC-MS analysis

3.2.1.1. LINEX system

- ♦OPTIC3 injector (ATAS GL)
- ♦FOCUS XYZ sample preparation robot (ATAS GL)

3.2.1.2. GC-MS system

- ♦Model 6890 GC / 5973 MSD (Agilent Technologies)
- ♦GC 353M / Varian 1200 (GL Sciences, Varian)
- ♦Trace GC / Poraris Q (Thermo Electron)

3.2.1.3. Column

- ♦Inert Cap Pesticides, 0.25mm I.D. x 30M (GL Sciences)

3.3. Analytical conditions

3.3.1. Injection (DMI: Difficult Matrix Introduction)

Initial temp.: 65°C, ramp rate: 5°C/sec, final temp.: 280°C, iso time: 5min, cool down: 65°C

3.3.2. Oven program

- ♦79°C (2min)-10°C/min-280°C (5min)




Figure 5. Optimized conditions in LINEX-DMI-GC/MS method




Figure 6. Repeated determinations of the normal hydrocarbons C10 to C30.

Table 1. Relative standard deviation of analytical results for the normal hydrocarbons

Hydrocarbon	C10	C12	C14	C16	C18	C20	C22	C24	C26	C28	C30
CV(%)	3.31	1.92	1.59	2.14	2.62	2.06	2.00	3.33	2.97	3.33	4.14

4. Evaluation

An internal standard method and peak-area ratio of the pesticide sample to the internal sample was used in the evaluation of the LINEX-DMI system. The ratio is also utilized as a compensation factor to calibrate analytical reproducibility.

$$[\text{Peak area ratio}] = \frac{[\text{Pesticide peak area}]}{[\text{Internal standard peak area}]}$$

5. Results and Discussion

5.1 Reproducibility of analytical results for the normal hydrocarbons

The determination of the normal hydrocarbons C10 to C30 was repeatedly carried out using the internal standard method. The reproducibility (n=45) was calibrated with the compensation factor, and the results are shown in Figure 6.

As the above results, the performance of PTV injector used in the study is satisfactory.

5.2 Repeated determinations of the pesticide without the LINEX-DMI

The standard pesticide samples diluted with acetone were repeatedly introduced into the capillary GC/MS system without any exchanging of the liner. The data for each compound were compensated with internal standard method.

As shown in Figure 7, the analytical values decreased as the repeated times increased. The reasons for this are that thermal decomposition and adsorption occurred inside the liner due to the use of acetone which can lead to an increase in the activity of the liner and decrease in the recovery even if using a de-active liner.

During the analysis, the compounds 2, 6 and 7 shown in Table 2 were partly decomposed to the compounds a, b and c, respectively. Therefore, endrin ketone and endrin ketone increased as the repeated times increased due to thermal decomposition of the original compound (endrin).




Figure 7. Repeated determinations of the pesticide

Table 2 Reproducibility (n=60) of the determination for standard sample 2

Pesticide	CV(%)
1 Fenobucarb	6.57
2 Bendiocarb	11.47
3 Dimethoate	11.96
4 Lindane	2.52
5 Diazinon	2.16
6 Ethiofenox	14.08
7 Carbaryl	14.35
8 Fenitrothion	9.92
9 Methiocab	16.73
10 Aldrin	4.56
11 Fenitron	5.58
12 Dieldrin	3.96
13 Endrin	14.81
14 p,p'-DDT	7.04
15 Ethofenox	7.77




Fig.8 Chromatograms for pesticides

Table 3. Decomposed products

Compounds
a Bendiocarb (decomposition)
b Ethiofenox (decomposition)
c Carbaryl (decomposition)
d Endrin ketone
e Endrin aldehyde

5.3 Repeated determination of pesticides with the LINEX-DMI

When the repeated determinations (25 times) were carried out, one liner was used for 5 runs and automatically exchanged with a new one by the LINEX.




Figure 9. Analytical results based on the LINEX

Table 3. Average peak-area ratio (cycle using a same liner) and relative standard deviation (%)

	1 cycle		2 cycle		3 cycle		4 cycle		5 cycle		Cycle average CV(%)
	AVG	CV(%)	AVG	CV(%)	AVG	CV(%)	AVG	CV(%)	AVG	CV(%)	
Endosulfan	0.74	0.84	0.70	2.07	0.71	2.97	0.70	1.13	0.73	3.07	2.84
Endosulfan S	0.54	1.47	0.47	4.03	0.52	4.90	0.48	1.50	0.53	3.96	6.19
Endosulfan R	0.43	2.15	0.39	2.97	0.41	4.29	0.39	1.06	0.41	4.83	4.73
Endosulfan Q	0.82	0.41	0.84	2.13	0.82	0.62	0.82	1.00	0.83	0.63	1.63
Endosulfan X	0.74	1.45	0.67	3.28	0.72	4.15	0.66	2.29	0.75	5.86	5.61
Endosulfan Y	0.44	2.77	0.43	3.09	0.41	8.52	0.35	2.31	0.42	8.01	8.35
Endosulfan Z	0.46	2.36	0.44	2.53	0.44	6.27	0.39	2.06	0.44	6.45	6.68
Endosulfan AA	0.46	1.72	0.40	3.91	0.43	5.55	0.40	1.48	0.44	6.75	5.78
Endosulfan AB	0.43	2.76	0.42	2.72	0.41	2.48	0.45	2.36	0.43	2.53	2.69
Aldrin	0.75	0.84	0.75	2.26	0.75	1.42	0.72	1.35	0.72	2.29	2.17
Endosulfan AC	0.71	1.28	0.66	2.48	0.70	2.09	0.66	1.32	0.69	3.05	3.05
Endosulfan AD	0.59	0.72	0.56	1.87	0.58	1.36	0.56	1.69	0.59	2.61	2.59
Endosulfan AE	0.56	2.83	0.54	4.58	0.54	6.64	0.56	1.93	0.56	5.34	2.48
Endosulfan AF	0.87	2.31	0.86	2.28	0.87	3.54	0.85	1.79	0.87	2.53	0.63
Endosulfan AG	0.77	4.58	0.74	4.54	0.73	7.95	0.75	3.38	0.77	5.18	6.82

As shown in Fig. 9, the recovery is poor for the first determination after changing the liner. This exhibits that it is necessary to condition a new liner. If removing the first data, the reproducibility was satisfactory for analytical results of each liner.

There are some considerable requirement and interest in throughput and matrix effect. The LINEX system is capable of making a variety of the liner-exchanging programs according to the requirements.

5.4 Preparation and determination based on the LINEX-DMI system

The preparation procedures of sesame are shown in Figure 10 (refer poster No. III-1-39C (SPW7-2) in detailed).

Sample 20 g
(Sesame)
 ↓ add water (20 mL) swelling at 30 min
Extract with acetone/n-hexane (1:2) 100 mL
 ↓ mix 30 sec by homogeniser
 ↓ and mix 20 min by reciprocate shaker
Centrifuge 5 min. at 1200 x g
 ↓
Decant supernatant
 ↓ → Residue was reextract with 50 mL
Evaporate
 ↓
Charge
 ↓ concentrated extract was
 ↓ diluted 20 mL with Hex
Elute
 ↓ eluted with 70 mL ACN
 ↓
Evaporate
 ↓ redissolve in 2 ml acetone
 ↓
GC/MS analysis

Figure 10. The preparation procedures of sesame

Exchanging of the liner was carried out for each run. 5µl extract obtained by the above preparation procedures and 0.3-ppm standard sample were introduced into the GC system by the LINEX-DMI, respectively. Figure 11 shows the chromatograms for the repeated determinations. Figure 12 shows a photograph of the liner in which sesame refuse remained after injection. The reproducibility of the determination of pesticides is reported in Table 5.




Figure 11. Determination of pesticides in sesame extract




Figure 12. Sesame refuse remained in the liner

Table 5. Reproducibility of the determinations for the pesticides in sesame extract

Pesticide	1 cycle	2 cycle	3 cycle	4 cycle	5 cycle	6 cycle	7 cycle	8 cycle	9 cycle	10 cycle	11 cycle	12 cycle	13 cycle	14 cycle	15 cycle	16 cycle	17 cycle	18 cycle	19 cycle	20 cycle	21 cycle	22 cycle	23 cycle	24 cycle	25 cycle																					
Endosulfan	0.74	0.84	0.70	2.07	0.71	2.97	0.70	1.13	0.73	3.07	2.84	0.74	0.84	0.70	2.07	0.71	2.97	0.70	1.13	0.73	3.07	2.84	0.74	0.84	0.70	2.07	0.71	2.97	0.70	1.13	0.73	3.07	2.84	0.74	0.84	0.70	2.07	0.71	2.97	0.70	1.13	0.73	3.07	2.84		
Endosulfan S	0.54	1.47	0.47	4.03	0.52	4.90	0.48	1.50	0.53	3.96	6.19	0.54	1.47	0.47	4.03	0.52	4.90	0.48	1.50	0.53	3.96	6.19	0.54	1.47	0.47	4.03	0.52	4.90	0.48	1.50	0.53	3.96	6.19	0.54	1.47	0.47	4.03	0.52	4.90	0.48	1.50	0.53	3.96	6.19		
Endosulfan R	0.43	2.15	0.39	2.97	0.41	4.29	0.39	1.06	0.41	4.83	4.73	0.43	2.15	0.39	2.97	0.41	4.29	0.39	1.06	0.41	4.83	4.73	0.43	2.15	0.39	2.97	0.41	4.29	0.39	1.06	0.41	4.83	4.73	0.43	2.15	0.39	2.97	0.41	4.29	0.39	1.06	0.41	4.83	4.73		
Endosulfan Q	0.82	0.41	0.84	2.13	0.82	0.62	0.82	1.00	0.83	0.63	1.63	0.82	0.41	0.84	2.13	0.82	0.62	0.82	1.00	0.83	0.63	1.63	1.63	0.82	0.41	0.84	2.13	0.82	0.62	0.82	1.00	0.83	0.63	1.63	1.63	0.82	0.41	0.84	2.13	0.82	0.62	0.82	1.00	0.83	0.63	1.63
Endosulfan X	0.74	1.45	0.67	3.28	0.72	4.15	0.66	2.29	0.75	5.86	5.61	0.74	1.45	0.67	3.28	0.72	4.15	0.66	2.29	0.75	5.86	5.61	0.74	1.45	0.67	3.28	0.72	4.15	0.66	2.29	0.75	5.86	5.61	0.74	1.45	0.67	3.28	0.72	4.15	0.66	2.29	0.75	5.86	5.61		
Endosulfan Y	0.44	2.77	0.43	3.09	0.41	8.52	0.35	2.31	0.42	8.01	8.35	0.44	2.77	0.43	3.09	0.41	8.52	0.35	2.31	0.42	8.01	8.35	0.44	2.77	0.43	3.09	0.41	8.52	0.35	2.31	0.42	8.01	8.35	0.44	2.77	0.43	3.09	0.41	8.52	0.35	2.31	0.42	8.01	8.35		
Endosulfan Z	0.46	2.36	0.44	2.53	0.44	6.27	0.39	2.06	0.44	6.45	6.68	0.46	2.36	0.44	2.53	0.44	6.27	0.39	2.06	0.44	6.45	6.68	0.46	2.36	0.44	2.53	0.44	6.27	0.39	2.06	0.44	6.45	6.68	0.46	2.36	0.44	2.53	0.44	6.27	0.39	2.06	0.44	6.45	6.68		
Endosulfan AA	0.46	1.72	0.40	3.91	0.43	5.55	0.40	1.48	0.44	6.75	5.78	0.46	1.72	0.40	3.91	0.43	5.55	0.40	1.48	0.44	6.75	5.78	0.46	1.72	0.40	3.91	0.43	5.55	0.40	1.48	0.44	6.75	5.78	0.46	1.72	0.40	3.91	0.43	5.55	0.40	1.48	0.44	6.75	5.78		
Endosulfan AB	0.43	2.76	0.42	2.72	0.41	2.48	0.45	2.36	0.43	2.53	2.69	0.43	2.76	0.42	2.72	0.41	2.48	0.45	2.36	0.43	2.53	2.69	2.69	0.43	2.76	0.42	2.72	0.41	2.48	0.45	2.36	0.43	2.53	2.69	2.69	0.43	2.76	0.42	2.72	0.41	2.48	0.45	2.36	0.43	2.53	2.69
Aldrin	0.75	0.84	0.75	2.26	0.75	1.42	0.72	1.35	0.72	2.29	2.17	0.75	0.84	0.75	2.26	0.75	1.42	0.72	1.35	0.72	2.29	2.17	0.75	0.84	0.75	2.26	0.75	1.42	0.72	1.35	0.72	2.29	2.17	0.75	0.84	0.75	2.26	0.75	1.42	0.72	1.35	0.72	2.29	2.17		
Endosulfan AC	0.71	1.28	0.66	2.48	0.70	2.09	0.66	1.32	0.69	3.05	3.05	0.71	1.28	0.66	2.48	0.70	2.09	0.66	1.32	0.69	3.05	3.05	0.71	1.28	0.66	2.48	0.70	2.09	0.66	1.32	0.69	3.05	3.05	0.71	1.28	0.66	2.48	0.70	2.09	0.66	1.32	0.69	3.05	3.05		
Endosulfan AD	0.59	0.72	0.56	1.87	0.58	1.36	0.56	1.69	0.59	2.61	2.59	0.59	0.72	0.56	1.87	0.58	1.36	0.56	1.69	0.59	2.61	2.59	0.59	0.72	0.56	1.87	0.58	1.36	0.56	1.69	0.59	2.61	2.59	0.59	0.72	0.56	1.87	0.58	1.36	0.56	1.69	0.59	2.61	2.59		
Endosulfan AE	0.56	2.83	0.54	4.58	0.54	6.64	0.56	1.93	0.56	5.34	2.48	0.56	2.83	0.54	4.58	0.54	6.64	0.56	1.93	0.56	5.34	2.48	2.48	0.56	2.83	0.54	4.58	0.54	6.64	0.56	1.93	0.56	5.34	2.48	2.48	0.56	2.83	0.54	4.58	0.54	6.64	0.56	1.93	0.56	5.34	2.48
Endosulfan AF	0.87	2.31	0.86	2.28	0.87	3.54	0.85	1.79	0.87	2.53	0.63	0.87	2.31	0.86	2.28	0.87	3.54	0.85	1.79	0.87	2.53	0.63	0.63	0.87	2.31	0.86	2.28	0.87	3.54	0.85	1.79	0.87	2.53	0.63	0.63	0.87	2.31	0.86	2.28	0.87	3.54	0.85	1.79	0.87	2.53	0.63
Endosulfan AG	0.77	4.58	0.74	4.54	0.73	7.95	0.75	3.38	0.77	5.18	6.82	0.77	4.58	0.74	4.54	0.73	7.95	0.75	3.38	0.77	5.18	6.82	6.82	0.77	4.58	0.74	4.54	0.73	7.95	0.75	3.38	0.77	5.18	6.82	6.82	0.77	4.58	0.74	4.54	0.73	7.95	0.75	3.38	0.77	5.18	6.82

6. Conclusion

The Contamination of the capillary column and the MS detector can be prevented by the application of the LINEX-DMI system.

The Ministry of Health, Labour and Welfare has established general provisions concerning the introduction of the "positive list" system for residual agricultural chemicals in foods in the year. For this, we can provide a methodology that is capable of achieving more effective and simpler preparation, higher throughput and reliability.