


Sample preparation for pesticides in fruits and vegetables

QuEChERS



EN 15662 Method

AOAC 2007.01 Method

GL Sciences - QuEChERS kit

QuEChERS (pronounced Catchers) is a sample preparation technique.

QuEChERS is an acronym for *Quick, Easy, Cheap, Effective, Rugged* and *Safe* which exactly explains its purpose. The kit makes the process of detecting pesticides in food samples easier and faster. It has become a common knowledge that pesticides influence the health of people. Different regulations all over the world make it important for the food industry to detect the different pesticides used on their import products. Therefore the QuEChERS kit was developed in 2003 and have become more important ever since.

There are 2 different methods written by both a European authority and an American authority referred to as the EN 15662 Method (EU) and the AOAC 2007.01 Method (USA).



How do they work?

The functionality of the QuEChERS kit is to extract any interfering matrices from the sample. The first step (explained on the next page) transfers the pesticides from the water to the added acetonitrile phase by applying MgSO_4 .

Lipids, sugars and pigments are the most important interfering matrices. The second step gets rid of these. Again MgSO_4 is added to separate any water. $\text{P}_{\text{primary}}\text{S}_{\text{secondary}}\text{A}_{\text{mine}}$ removes fatty acids and sugars, $\text{G}_{\text{raphite}}\text{C}_{\text{arbon}}\text{B}_{\text{lack}}$ removes pigments and sterols. C18 makes sure that the pesticides do not dissolve in the lipid layer that arises in fat fruits and vegetables (sample > 1% fat).

The extraction can be done after homogenizing 10 grams of the sample. When working with dried fruits, add some water before homogenizing. Apply the sample to the provided 50 mL test tubes for the first step. Starting the second step can be done by applying the aliquot supernatant to the provided 2 mL or 15 mL centrifuge tubes.


What QuEChERS kit to use?

On the next page one can find the Buffered Extraction Kits. Pay attention to choosing the correct one for the international market you're producing for (either EU or USA).

On the back of this folder a list of QuEChERS kits can be found with a description on what fruits and vegetables they should be used with. Again observe the EU and USA method being displayed.

Using the wrong QuEChERS kit can lead to a contaminated chromatogram as well as a contaminated column. Make sure you do know what matrices inside the fruits you are dealing with.

Method, Step 1



GL Sciences B.V. Buffered Extraction Kit

EN 15562 Method Buffered Extraction Kit

4 g MgSO₄, 1 g NaCl, 1 g NaCitrate, 0.5 g diSodium Citrate Sesquihydrate


P/No: 8010-1001

AOAC 2007.01 Method Buffered Extraction Kit

6 g MgSO₄, 1.5 g NaCitrate

P/No: 8010-1002

Method, Step 2



GL Sciences B.V. can supply you with the best columns for both GC-MS and LC-MS.

InertCap® GC Columns



InertSustain® LC Columns



Inertsil® LC Columns



Order Information

	Matter Removed from sample	Examples	EN 15562 Method	AOAC 2007.01 METHOD
	General Purpose - Polar Organic Acids - Some Sugars - Some Lipids	• Apple • Peach • Broccoli • Fresh Peas	25 mg PSA, 150 mg MgSO ₄ 2 ml - 100/box P/No: 8010-1003	50 mg PSA, 150 mg MgSO ₄ 2 mL - 100/box P/No: 8010-1005
	Fruits and vegetables with fats and waxes - Polar Organic Acids - Some Sugars - Lipids - Sterols	• Oranges • Lemons • Animal Sourced • Cereals	25 mg PSA, 25 mg C18 150 mg MgSO ₄ 2 mL - 100/box P/No: 8010-1007	50 mg PSA, 50 mg C18 150 mg MgSO ₄ 2 mL - 100/box P/No: 8010-1009
	Pigmented fruits and vegetables* - Polar Organic Acids - Some Sugars - Some Lipids - Carotinoïdes and Chlorophyll	• Zucchini • Strawberry • Dried Peas • Tomatoes	25 mg PSA, 2.5 mg GCB 150 mg MgSO ₄ 2 mL - 100/box P/No: 8010-1011	50 mg PSA, 50 mg GCB 150 mg MgSO ₄ 2 mL - 100/box P/No: 8010-1013
	Highly pigmented fruits and vegetables* - Polar Organic Acids - Some Sugars - Some Lipids - Carotinoïdes and Chlorophyll	• Black Berries • Spinach • Red Peppers	25 mg PSA, 7.5 mg GCB 150 mg MgSO ₄ 2 mL - 100/box P/No: 8010-1015	
	Fruits and vegetables with pigments and fats* - Polar Organic Acids - Some Sugars - Some Lipids - Carotinoïdes and Chlorophyll	• Dairy products • Curry Spices • Baby food		50 mg PSA, 50 mg C18 50 mg GCB, 150 mg MgSO ₄ 2 mL - 100/box P/No: 8010-1017
	Other food methods - Fats - Lipids - Proteins	• Olive Oil • Citrus Oil • Liver		400 mg PSA, 400 mg C18 400 mg GCB, 1200 mg MgSO ₄ 15 mL - 50/box P/No: 8010-1018
	All food types - All matrix interfering materials			25 mg C18 150 mg MgSO ₄ 2 mL - 100/box P/No: 8010-1019
				150 mg C18 900 mg MgSO ₄ 15 mL - 50/box P/No: 8010-1020
				50 mg PSA, 50 mg C18 7.5 mg GCB, 150 mg MgSO ₄ 2 mL - 100/box P/No: 8010-1021
				400 mg PSA, 400 mg C18 45 mg GCB, 1200 mg MgSO ₄ 15 mL - 50/box P/No: 8010-1022

* Not for use with planar pesticides

Contact

GL Sciences B.V.
De Sleutel 9
5652AS Eindhoven
The Netherlands
Tel. +31 (0)40 2549531
info@glsciences.eu
www.glsciences.eu